DENTSPLY

CERAMCO®

*Directions For Use
Mode d'emploi
Istruzioni per l'uso
Gebrauchsanleitung
Instrucciones de uso*

¹Vita® and 3D-Master® are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co. ²Data on file.
¹Vita® et 3D-Master® sont des marques déposées de Vita Zahnfabrik H. Rauter GmbH & Co. ²Données sur fichier.
¹Vita® e 3D-Master® sono marchi registrati di Vita Zahnfabrik H. Rauter GmbH & Co. ²Dati registrati.
¹Vita® und 3D-Master® sind eingetragene Warenzeichen der Vita Zahnfabrik H. Rauter GmbH & Co. ²Daten im Archiv.
¹Vita® y 3D-Master® son marcas registradas de Vita Zahnfabrik H. Rauter GmbH & Co. ²Datos en archivo.

Manufacturer:
DENTSPLY International Inc.
DENTSPLY Ceramco
570 West College Ave.
York, PA 17401 U.S.A.
1-800-243-1942
prosthetics.dentsply.com

ECIREP DeguDent GmbH
Rodenbacher Chaussee 4
63457 Hanau-Wolfgang
Germany
Telefon +49/6181/5950

Distributed by:
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3
905-968-5374

A
Academy
of General Dentistry
PROUD SPONSOR

CE
0086
Rx only

Ceramco3 is a complete porcelain-fused-to-metal system that has been engineered to provide exceptional handling, thermal stability and esthetic properties. The entire system, from the opaques through the dentins, natural enamels and opal enamels impart the fluorescence of natural dentition.

INDICATIONS

The Ceramco3 veneering porcelains may be used for all single and multiple unit porcelain fused-to-metal fixed prosthodontics using conventional high-fusing ceramic alloys with or without silver.

CONTRAINDICATIONS

Only the indications listed above are suitable.

WARNINGS/ADVERSE REACTIONS

Consult Ceramco3 SDS (Safety Data Sheets) for safe handling and usage. Use only with recommended Ceramco 3 system components. Use protective equipment to guard against dust when grinding.

PRECAUTIONS

Refer to the special notes throughout the instructions for best practices for using Ceramco3 materials.

INSTRUCTIONS

PORCELAIN FURNACE

For optimum results, be sure your porcelain furnace is accurately calibrated and manufacturer recommended firing times and temperatures are followed. If necessary, adjust the furnace parameters to obtain the correct visual indicator.

SAGGER TRAYS

During the firing cycle, large sagger trays will absorb some of the heat and may effect the appearance of the porcelain. For optimum results, only use the honeycomb style trays with the metal firing pins.

LIQUIDS

Caution should be used in selecting modeling and staining liquids. DENTSPLY liquids are recommended for their ability to improve the handling and carving of the porcelain, to reduce slumping during build-up, and to prevent any silver discoloration. Distilled water may also be used. However, other build-up and/or stain liquids may leave a residue behind that could cause discoloration.

STORAGE

Store materials in packages with the lid tightly sealed. Do not leave open bottles where dust and debris may contaminate the porcelain. Replace cap when not in use. Avoid vibration and heat. Keep dry. Shake well before use. The lot number is printed on each bottle; please refer to the lot number in all correspondence. Do not use after expiration date. DO NOT FREEZE, Margin & Opaque Correction Liquid, REF # 301916.

ALLOY SELECTION

Ceramco3 is compatible with high noble, noble and predominately base alloys. Be sure to consult the alloy manufacturer for alloy compositions and coefficient of thermal expansion data. Ceramco3 will resist discoloration on silver containing alloys; however, routine purging of the porcelain furnace and sagger trays is highly recommended.

(This refers to the cooling cycle programmable into most porcelain furnaces.)

CTE values listed above in units of $\times 10^{-6} \text{K}^{-1}$ (20-600 °C)

The Ceramco3 porcelain system includes a paste opaque and/or a powder opaque in shades keyed to Vita[®] classical and Vita[®] 3D-Master[®] Shade Series², as well as DENTSPLY's Illuminé[™] shades. Refer to the specific section in this manual for application and firing guidelines.

PASTE OPAQUE APPLICATION

1. Oxidize or degass and clean the alloy according to the alloy manufacturer's instructions.
2. Rinse a flat brush in distilled water and wipe it on a dry paper towel. Check all castings to make sure that all surfaces are dry.
3. Select the proper base paste depending on the shade. Regular base paste is used for the Vita shades and white base paste is used for the Illuminé shades.
4. Dispense a small amount of the base paste onto a pallet. Water must not be mixed into the paste. If necessary, the paste may be thinned with Opaque Modifier Fluid. (DO NOT thin the Illuminé base paste.) Use a brush and a light touch to apply as even a paste coat as possible onto the metal substructure. Achieve an even coat and do not allow the material to build up or puddle. For the Illuminé base paste, vibrate to allow the base paste to smooth out. Brush strokes that are present will flatten out with light vibration from a serrated instrument.

Base Paste Application

Note: At this point the metal should be masked out. An application of the base paste that is too thin will result in shades that are too gray or low in value.

5. Use a hemostat and hold the restoration over the appropriate collecting bowl. Using the following chart, sprinkle the correct crystals onto the base paste.

Shades		Base Paste	Crystals used After Base Paste	Crystals used After Shaded Paste
i1, i2, i3, i4, i5, i6, i7, i8		Illuminé White	White	White
A1, A2, A3 A3.5, A4		Regular Regular	Medium Dark	Light Light
B1, B2 B3, B4		Regular Regular	Medium Dark	Light Light
C1 C2, C3, C4		Regular Regular	Medium Dark	Light Light
D2, D4 D3		Regular Regular	Medium Dark	Light Light
Ceramco®3 C3 Shades Series	Keyed to Vita®3D Master®⁰¹ Shades²	Base Paste	Crystals used After Base Paste	Crystals used After Shaded Paste
C3-101	1M1	Illuminé White	Medium	Light
C3-102	1M2	Illuminé White	Medium	Light
C3-203	2L1.5	Regular	Medium	Light
C3-204	2L2.5	Regular	Medium	Light
C3-205	2M1	Regular	Medium	Light
C3-206	2M2	Regular	Medium	Light
C3-207	2M3	Regular	Medium	Light
C3-208	2R1.5	Regular	Medium	Light
C3-209	2R2.5	Regular	Medium	Light
C3-310	3L1.5	Regular	Dark	Light
C3-311	3L2.5	Regular	Dark	Light
C3-312	3M1	Regular	Medium	Light
C3-313	3M2	Regular	Medium	Light
C3-314	3M3	Regular	Dark	Light
C3-315	3R1.5	Regular	Dark	Light
C3-316	3R2.5	Regular	Dark	Light
C3-417	4L1.5	Regular	Dark	Light
C3-418	4L2.5	Regular	Dark	Light
C3-419	4M1	Regular	Dark	Light
C3-420	4M2	Regular	Dark	Light
C3-421	4M3	Regular	Dark	Light
C3-422	4R1.5	Regular	Dark	Light
C3-423	4R2.5	Regular	Dark	Light
C3-524	5M1	Regular	Dark	Light
C3-525	5M2	Regular	Dark	Light
C3-526	5M3	Regular	Dark	Light

1. Vita® and 3D-Master® are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co. 2. Data on file.

Crystal Application

6. Tap the hemostat very lightly to remove excess crystals. Do Not Vibrate! Use a dry brush to remove any crystals from the inside of the coping and any non-porcelain bearing areas. Use a brush at a 45° angle to remove any crystals that may remain in the margin area.
7. Dry and fire according to the recommended temperatures. The correct visual indicator is a sandpaper-like finish.

First Firing

Note: Drying too rapidly results in small voids in the base paste under the crystals. These small voids may come to the surface during subsequent firings. Following the recommended dry and pre-dry times will ensure this does not occur. It is also effective to place the restoration on a sagger tray and put the assembly on a warm surface or hot plate prior to placing in the porcelain furnace.

8. Dispense a small amount of the prescribed shaded paste on a pallet. Using moderate pressure paint the base paste with an even layer of the shaded paste. A light vibration or tapping should be used to smooth out the material.

Shaded Paste Application

9. If paste opaque modifiers are required, apply at this time. They may be used in full concentration or mixed with the shaded paste. Remove excess paste from the inside of the coping.

Modifier Application

10. Hold the restoration with a hemostat over the collecting bowl and dust all painted surfaces with the light or white crystals depending on the shade. Refer to the table on page 3. Tap the hemostat lightly to remove excess crystals from the coping.

Crystals after the Shaded Paste

11. To achieve an extremely thin, yet opaque layer at the marginal areas, use a soft, dry, sable brush to draw off some of the crystals by brushing at a 45° angle to the margin. Do not use heavy pressure.
12. Dry and fire according to the recommended temperatures. The correct visual indicator is a sandpaper-like finish.

Second Firing

POWDER OPAQUE APPLICATION

First Opaque Porcelain Application

1. Oxidize or degass and clean the alloy according to the alloy manufacturer's instructions.
2. Mix the opaque powder with DENTSPLY Build-up Modeling Liquid E or U, or distilled water to a creamy consistency.
3. Apply the first opaque porcelain layer thinly and evenly with a brush or instrument.

First Opaque Application

4. Tap the hemostat very lightly to condense and smooth the opaque porcelain surface. Do not use excessive vibration or allow thick layers of opaque porcelain to accumulate in occlusal, interproximal or marginal areas. Thick layers of opaque porcelain will tear or cause a fissure during the firing process.
5. Dry the opaque porcelain with warm, forced-air equipment or by holding the casework close to, but not in, the entrance of the furnace. Do not use extreme heat to dry the opaque layer. Excessive heat causes blisters, peeling or lifting of the opaque porcelain.
6. Dry and fire according to the recommended temperatures. Correct visual indicator: Slight sheen.

First Opaque Fired

Second Opaque Porcelain Application

1. Moisten the dried opaque porcelain mixture with distilled water only. The consistency for the second layer of opaque porcelain should be slightly thicker than the initial mix.
2. Apply the second opaque porcelain layer using the same technique as for the first layer. Make sure this layer is thick enough to completely mask the metal substructure.

Second Opaque Application

3. If opaque modifiers are required, apply them at this time.
4. Do not allow thick layers of opaque porcelains to accumulate on concave areas such as the occlusal grooves, interproximal areas or lingual shoulders. Excessive thickness of the opaque layer will fissure during firing.
5. Dry and fire according to the recommended temperatures. Correct visual indicator: Slight sheen.

Second Opaque Firing

The porcelain margin can be an integral part of a premium restoration, maximizing esthetics and oral health. The porcelain margin helps avoid metal visibility, tissue discoloration and possible overbuilding at the gingival. A patented carrier medium developed for Ceramco3 margin porcelain facilitates esthetic, easy, accurate and rapid placement for predictable - and profitable porcelain margins.

GENERAL PRINCIPLES

1. Dispense only enough material for immediate use to avoid waste. Dried Ceramco3 margin porcelain prepared with the DENTSPLY Margin and Opaque Correction Liquid can be reused. Select the correct margin porcelain from the following chart. Mixing some of the margin modifier porcelain into the shaded margin porcelain can easily change the porcelain shades.

Margin Porcelain	Corresponding Shade	Ceramco [®] 3 C3 Shades Series	Keyed to Vita [®] 3D Master [®] Shades ²
i1	i1, i3		
i2	i2, i4		
i5	i5, i6, i7, i8		
A1	A1, A2	C3-101	1M1
		C3-102	1M2
		C3-206	2M2
		C3-209	2R2.5
A3	A3, A3.5 B3 D3, D4	C3-204	2L2.5
		C3-207	2M3
		C3-311	3L2.5
		C3-313	3M2
		C3-314	3M3
		C3-315	3R1.5
		C3-316	3R2.5
		C3-419	4M1
		C3-422	4R1.5

Margin Porcelain	Corresponding Shade	Ceramco [®] 3 C3 Shades Series	Keyed to Vita [®] 3D Master [®] Shades ²
A4	A4 B4	C3-418	4L2.5
		C3-420	4M2
		C3-421	4M3
		C3-423	4R2.5
		C3-525	5M2
B2	B1, B2	C3-526	5M3
		C3-203	2L1.5
C1	C1, C2 D2	C3-205	2M1
		C3-208	2R1.5
		C3-312	3M1
C3	C3, C4	C3-310	3L1.5
		C3-417	4L1.5
		C3-524	5M1

2. Clean brushes with water after each margin application to prevent contamination.
3. Complete the 5 minute dry and 5 minute preheat cycles to remove organic material and avoid discoloration. For furnaces without a "preheat", dry for 10 minutes.
4. Fire the porcelain on honeycomb sagger trays with metal pins.
5. Properly fired Ceramco3 margin porcelain will have a shiny yet grainy appearance. Firing the dentin porcelain higher than the margin porcelain will result in rounded margins.
6. Be sure the fired porcelain has the correct visual indicator.

TOOTH MARGIN PREPARATION

1. A 90° shoulder preparation or deep chamfer is recommended.
2. Normal tooth reduction (1.5mm) is recommended.

SUBSTRUCTURE DESIGN

Extend the facial surface of the metal substructure so that the metal ends at the axial gingival corner of the margin preparation or extends slightly onto the gingival ledge.

Proper Metal Design

OPAQUE PROCEDURE

Apply and fire the Ceramco3 paste opaque or Ceramco3 powder opaque porcelain using your normal procedures.

Be sure the opaque completely covers the metal at the gingival margin. Failure to do so will make the completed porcelain margin appear gray.

DIE SEALER APPLICATION

1. Using DENTSPLY Prosthetics die sealer, apply a thin, even layer of sealer over the die margin area. Blow off any excess die sealer and allow the liquid to thoroughly dry.
2. Reapply the die sealer at least 2 more times and allow to dry.
3. Mark the facial margin with a wax pencil.

Die Sealer Application

DIE RELEASE APPLICATION

1. Be sure the die margin is clean.
2. Using DENTSPLY Prosthetics die release, apply a thin even layer of die release over the margin area. Blow off the excess die release.
3. Apply and allow to dry several more times until the die has a very slight sheen.

Die Release Application

FIRST MARGIN APPLICATION

1. Place a small amount of the shaded margin porcelain on a clean palette.
2. If modifier is necessary, add the modifier porcelain to the shaded margin porcelain.
3. Mix the Ceramco3 Margin and Opaque Correction Liquid with the Margin Porcelain to a thick creamy consistency. This mixture will air harden in approximately 5 minutes. If additional working time is required, mix in a small amount of distilled water.
4. Place the porcelain mixture in the gingival area, pushing the mixture down to but not over the gingival margin. Do not over build.

First Shoulder Application

Note: All porcelains shrink by volume. The larger the porcelain mass the more apparent the firing shrinkage. Better results are achieved by applying a small amount of porcelain on the first porcelain application.

5. Do not allow the restoration to become completely dry before lifting off the model.
6. Carefully lift the coping from the die.

First Shoulder Lift-off

Note: If inspection reveals the need for adjustment, carefully return the coping to the die and add fresh material. It is not necessary to remove the unfired porcelain.

7. Dry and fire the restoration according to the recommended temperatures.
8. The correct visual indicator will be shiny yet grainy.

First Shoulder Fired

SECOND MARGIN APPLICATION

1. Be sure the die is clean and reapply the die release. Allow to dry.
2. Place the coping on the die. Make a new mix of Ceramco3 Margin porcelain and Ceramco3 Margin and Opaque Correction Liquid. Apply the mixture to the gingival margin.

Second Shoulder Application

3. Do not allow the restoration to become completely dry before lifting off the model.

Second Shoulder Lift-Off

4. Dry and fire the restoration according to the recommended temperatures.
5. The correct visual indicator will be shiny yet grainy.

Shoulder Margin Completed

Note: If additional applications are required, repeat the procedures used for the second margin application.

Natural Enamel, Opal Enamel, and Mamelon Porcelain Application

GENERAL PRINCIPLES

1. Mix the Ceramco3 porcelains with DENTSPLY Build-up Modeling Liquids E or U, or distilled water to a paste-like consistency. Build-up liquids E or U will make the porcelain more "plastic-like" as compared to using distilled water. Caution must be used when selecting other build-up liquids. Some liquids may not completely burnout during the firing cycle and cause the porcelain to discolor or turn yellow/green on silver containing alloys.
2. Ceramco3 is a fine-grained porcelain. The degree to which the porcelain is condensed is based on individual techniques. Additional condensing removes liquid and packs the porcelain particles more thoroughly. This results in porcelain with a higher chroma and less shrinkage. Less condensing results in additional shrinkage and shades with lower chroma. Normally only slight condensing of the porcelain is necessary.
3. Apply DENTSPLY Prosthetics die sealer to the areas of the model that will come in contact with the porcelain. This will prevent the die stone from absorbing liquid from the porcelain.
4. Fire the porcelain on honeycomb sagger trays with metal pins.
5. Dry and fire according to the recommended temperatures.
6. Be sure the fired porcelain has the correct visual indicator.

CERAMCO3 OPACEOUS DENTINS


Natural teeth exhibit numerous color zones. These areas have various degrees of translucency, chroma and opacity. The experienced technician learns to use a variety of shades and modifiers to effectively mimic these characteristics. Frequently, the task is compromised by circumstances that do not allow for adequate porcelain thickness. The Ceramco3 Opaceous Dentin porcelains address the problem of limited space, simplify internal color modifications and aid in the overall build-up.

APPLICATIONS

The Ceramco3 Opaque Dentin porcelains are the same hue as the Ceramco3 dentin shades, but have approximately 10% more opacity. They may be used undiluted, or mixed with the Ceramco3 dentins or dentin modifiers.


Thin Areas

Frequently insufficient dentin/incisal thickness exists (less than 1.0mm) over the opaque porcelain. The shade in this thin area can be improved by applying a layer of the opaque dentin porcelain under the same dentin porcelain shade. When contouring, you will be able to reduce the area where the opaque dentin was layered and the prescribed shade will be maintained. Complete the building and firing of the restoration in the normal manner.


Facial Margin (Gingival Color)

If the metal substructure has been designed with a feather margin, the opaque porcelain around the facial margin leaves minimal room for the dentin porcelain, and the opaque porcelain will be highlighted. Esthetics will be improved by applying opaque dentin porcelain alone or as a blend with 10% to 20% of one of the opaque dentin modifiers. This technique creates a chromatic cervical third, common in natural dentition.


Pontics

A color imbalance frequently occurs between the pontics and the abutments of a multiple-unit fixed partial denture. The color difference is caused by the absence of a full metal substructure in the pontic area. Ceramco3 opaque dentin porcelain applied to the bottom and the cervical area of the pontic creates a uniform color across the fixed partial denture.

Lingual Fossas and Occlusal Surfaces

Often, inadequate space is available on the lingual surface of an anterior unit, or on the occlusal surface of a posterior unit. Ceramco3 opaque dentin porcelain applied alone or as a blend with 10% to 20% opaque dentin modifiers in these thin areas prevents the opaque porcelain from being highlighted.

CERAMCO3 DENTIN AND DENTIN MODIFIER PORCELAIN

The Ceramco3 Dentin Modifiers have the same hue and opacity as the Ceramco3 dentin porcelain but are approximately 70% higher in chroma (color saturation). The dentin modifiers can be used undiluted or mixed with the Ceramco3 Dentins, Opaceous Dentins, Natural Enamels, or Opal Enamels for a variety of effects. There are 6 New Shaded Dentin Modifiers and 6 New Hue Modifiers available. Select the recommended dentin modifier that corresponds with the desired shade according to the following chart. For a very subtle effect, mix in a very small amount of dentin modifier with the dentin porcelain.

New Ceramco 3 High Chroma (HC1-6) Dentin Modifiers

Mixing Scheme

Shade	Ratio
A1	100% HC2
A2	75% HC2, 25% HC5
A3	50% HC2, 50% HC5
A3.5	100% HC5
A4	75% HC5, 25% HC3
B1	100% HC1
B2	75% HC1, 25% HC4
B3	50% HC1, 50% HC4
B4	80% HC4, 20% HC3
C1	100% HC3
C2	75% HC3, 25% HC6
C3	50% HC3, 50% HC6
C4	100% HC6
D2	50% HC5, 50% HC3
D3	75% HC2, 25% HC4
D4	100% HC3

Mixing Scheme for Shade Series

Shade	Keyed to VITA® 3D-Master® Shade Series¹	Ratio
101	1M1	100% HC1
102	1M2	100% HC2
203	2L1.5	100% HC3
204	2L2.5	75% HC1, 25% HC4
205	2M1	50% HC2, 50% HC3
206	2M2	75% HC2, 25% HC5
207	2M3	50% HC2, 50% HC5
208	2R1.5	75% HC2, 25% HC5
209	2R2.5	75% HC2, 25% HC5
310	3L1.5	75% HC2, 25% HC4
311	3L2.5	50% HC1, 50% HC4
312	3M1	100% HC3
313	3M2	50% HC2, 50% HC5
314	3M3	80% HC4, 20% HC3
315	3R1.5	75% HC3, 25% HC6
316	3R2.5	50% HC2, 50% HC5
417	4L1.5	50% HC3, 50% HC6
418	4L2.5	80% HC4, 20% HC3
419	4M1	50% HC3, 50% HC6
420	4M2	100% HC5
421	4M3	80% HC4, 20% HC3
422	4R1.5	75% HC3, 25% HC6
423	4R2.5	75% HC5, 25% HC3
524	5M1	100% HC6
525	5M2	100% HC6
526	5M3	75% HC5, 25% HC3

CERAMCO3 NATURAL ENAMELS, OPAL ENAMELS AND MAMELON PORCELAIN

The Ceramco3 Natural Enamels are non-opalescent and can be used in thin or thick layers over the dentin to duplicate the vitality of natural teeth. They can also be mixed with the dentin or opal enamel porcelains for a variety of effects. The Ceramco3 SoftWear Enamels can be used in place of the Natural Enamels, to provide low wear characteristics against opposing natural dentition. The Ceramco3 Opal Enamels are an intense opal enamel porcelain. They can be mixed in with the Natural Enamels or used as an effect powder. The mamelon porcelains can be used in full concentration or diluted with the clear enamel porcelain. They are normally placed on the dentin porcelain lobes prior to the enamel porcelain application.

DENTIN BUILD-UP TECHNIQUE:

1. Areas on the model that the porcelain will contact should be completely sealed with DENTSPLY Prosthetics Die Sealer. Other sealers or oils may leave a residue that will cause discoloration and should not be used.
2. Mix the porcelains with DENTSPLY Build-up Modeling Liquids E or U, or distilled water to a paste-like consistency.

1. Vita® and 3D-Master® are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co. 2. Data on file.

3. Build dentin porcelain or opaque dentin/dentin modifier mixture around the cervical one third. Condense slightly. If building a multiple-unit fixed-partial denture, apply the dentin porcelain mixture to the saddle area of the pontic and seat the framework on the model. Continue to apply the dentin porcelain until full contour has been achieved. Condense slightly.

Opaque Dentin Application

4. Cutback the incisal one third to allow room for the enamel porcelains. If required, cutback to accentuate the lobe formation found in natural teeth. Be sure to remove adequate porcelain in the interproximal area.

Dentin Cut-Back

5. If mamelons are required, place a small amount of mamelon porcelain on the dentin lobes. To simulate natural dentition, vary the length of the mamelons.

Mamelon Application

6. Select the appropriate Natural Enamel Porcelain from Table 1. Apply the enamel porcelain to the incisal or occlusal one-third, overbuilding by approximately 10%. If desired layer the selected enamel porcelains in narrow vertical columns across the facial-incisal area to simulate the enamel rod formations found in natural teeth. Use Natural Enamel porcelains of different colors, the clear enamel porcelain, or the Opal Enamel porcelains arranged in discrete parallel columns to create contrast. Condense lightly.

6a. Table 1 can be used as a guide line for selecting an enamel porcelain. For the Illuminé shades, enamels will have a dramatic effect in the shade of the final restoration. As a rule of thumb:

1. Any translucent porcelain (Clear, Opal Clear) will make the shade lower in value (gray).
2. Due to the opalescence, any opal enamel will give the shade a very slight "warm" (yellow/orange) appearance.
3. To achieve a white/bright enamel, use Natural White, or Extra Light enamel.
4. To achieve a cold blue appearance, use Natural Blue enamel.

TABLE 1
SUGGESTED ENAMELS and MAMELONS

Shade	Natural Enamel	Opal Effect Enamel	Mamelon
i1	Extra Light	Opal Light	Yellow-Orange
i2	Extra Light	Opal Light	Yellow-Orange
i3	Extra Light	Opal Light	Yellow-Orange
i4	White	Opal White	Yellow-Orange
i5	White	Opal White	Yellow-Orange
i6	White	Opal White	Yellow-Orange
i7	White	Opal White	Yellow-Orange
i8	White	Opal White	Yellow-Orange
A1	Extra Light	Opal Light	Yellow-Orange
A2	Extra Light	Opal Light	Yellow-Orange
A3	Light	Opal Light	Red-Orange
A3.5	Light	Opal Medium	Red-Orange
A4	Light	Opal Medium	Red-Orange
B1	Extra Light	Opal White	Yellow-Orange
B2	Extra Light	Opal Light	Yellow-Orange
B3	Light	Opal Light	Red-Orange
B4	Light	Opal Light	Red-Orange
C1	Light	Opal Light	Yellow-Orange
C2	Light	Opal Light	Yellow-Orange
C3	Medium	Opal Medium	Red-Orange
C4	Dark	Opal Medium	Red-Orange
D2	Extra Light	Opal Light	Yellow-Orange
D3	Light	Opal Light	Yellow-Orange
D4	Medium	Opal Light	Red-Orange

Ceramco®3 C3 Shades Series	Keyed to Vita®3D Master® Shades ²	Enamel
C3-101	1M1	C-3-Light (EN-2)
C3-102	1M2	C-3-Light (EN-2)
C3-206	2M2	C-3-Light (EN-2)
C3-209	2R2.5	C-3-Light (EN-2)
C3-204	2L2.5	C-3-Light (EN-2)
C3-207	2M3	C-3-Light (EN-2)
C3-311	3L2.5	C-3-Light (EN-2)
C3-313	3M2	C-3-Light (EN-2)
C3-314	3M3	C-3-Light (EN-2)
C3-315	3R1.5	C-3-Light (EN-2)
C3-316	3R2.5	C-3-Light (EN-2)
C3-419	4M1	C-3-Light (EN-2)
C3-422	4R1.5	C-3-Light (EN-2)

Ceramco®3 C3 Shades Series	Keyed to Vita®3D Master® Shades ²	Enamel
C3-418	4L2.5	C-3-Light (EN-2)
C3-420	4M2	C-3-Light (EN-2)
C3-421	4M3	C-3-Light (EN-2)
C3-423	4R2.5	C-3-Light (EN-2)
C3-525	5M2	C-3-Light (EN-2)
C3-526	5M3	C-3-Light (EN-2)
C3-203	2L1.5	C-3-Light (EN-2)
C3-205	2M1	C-3-Light (EN-2)
C3-208	2R1.5	C-3-Light (EN-2)
C3-312	3M1	C-3-Light (EN-2)
C3-310	3L1.5	C-3-Light (EN-2)
C3-417	4L1.5	C-3-Light (EN-2)
C3-524	5M1	C-3-Light (EN-2)

1. Vita® and 3D-Master® are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co. 2. Data on file.

Basic Enamel Application

Clear Enamel overlay

Vertical Columns of Clear, Opal and Natural Enamels

Completed Build-Up

7. Remove the coping from the model. Apply a small amount of enamel porcelain to the mesial/distal contacts. Slightly condense the porcelain.
8. If building a multiple-unit fixed-partial denture, use a thin bladed instrument and cut slightly into the interproximals. It is not necessary to cut all the way down to the opaque.
9. Dry the build-up away from direct heat for at least five minutes before bringing it close to the open furnace entrance. Increase the drying time to ten minutes for large multiple-unit fixed-partial dentures.
10. Fire according to the recommended temperatures. A proper visual indicator is shiny with a small amount of surface texture.

Correct Visual Indicator for First Firing

CONTOURING & CLEANING

1. After firing, use fine-grit, non-contaminating stones, discs, or diamond burs to refine the anatomy. If no porcelain additions are needed, the case may be glazed after thorough cleaning.
2. If additions are required or before glazing, lightly blast the porcelain surface with unrecycled 50-micron aluminum oxide at 20 psi. Clean the case with a steam cleaner or ultrasonic bath using distilled water.

SECOND APPLICATION OF DENTIN & ENAMEL PORCELAINS

Apply the additional porcelain using the same steps and techniques described for the initial application. Match the layering of the opal enamel porcelain over the dentin porcelain. Second applications of Ceramco3 porcelains require the same care in drying and firing as the first layer.

OPAQUE CORRECTION

If during the contouring procedures the porcelain surface is removed and the metal surface is exposed, this may be easily repaired by using the opaque correction material. This correction opaque should be used for small minor repairs only.

1. Lightly blast the surface that is to be repaired with aluminum oxide at 20 psi. Then steam clean or clean in an ultrasonic cleaner in distilled water for 5 minutes.
2. Place a small amount either the light, medium or dark correction opaque onto a pallet. Mix with a small amount of Ceramco3 Margin/Opaque Correction Liquid to a paste like consistency.

3. Apply the opaque mixture to the exposed metal and allow to dry for about 5 minutes. The opaque will appear chalky and hard.
4. Apply either Ceramco3 opaque dentin, dentin or add-on porcelain and fire according to the recommended firing temperatures.

ADD-ON PORCELAIN

The Ceramco3 Add-On porcelain is available in a translucent enamel, dentin light, white, medium and dark, and tissue colored pink, reddish-pink, salmon, and dark. They fire at slightly lower temperatures than the regular dentin porcelain. These porcelains are always fired under vacuum and are used when a low temperature add-on is required.

1. Lightly blast the surface that is to be repaired with aluminum oxide at 20 psi. Then steam clean or clean in an ultrasonic cleaner in distilled water for 5 minutes.
2. Mix the porcelains with DENTSPLY Build-up Liquids E or U, or distilled water to a paste-like consistency.
3. Apply the mixture to the desired area.
4. Fire according to the recommended firing temperatures.

DENTSPLY

CERAMCO[®]3

Final Margin

FINAL MARGIN APPLICATION

The Final Margin is designed to be a final porcelain margin application that can be used to repair small margin imperfections. It is designed to be fired at the same time as the glaze cycle.

Final Margin Porcelain	Corresponding Shade	Ceramco [®] 3 C3 Shades Series	Keyed to Vita [®] 3D Master [®] Shades ²
i1	i1, i3		
i2	i2, i4		
i5	i5, i6, i7, i8		
A1	A1, A2	C3-101	1M1
		C3-102	1M2
		C3-206	2M2
		C3-209	2R2.5
A3	A3, A3.5 B3 D3, D4	C3-204	2L2.5
		C3-207	2M3
		C3-311	3L2.5
		C3-313	3M2
		C3-314	3M3
		C3-315	3R1.5
		C3-316	3R2.5
		C3-419	4M1
		C3-422	4R1.5

Final Margin Porcelain	Corresponding Shade	Ceramco [®] 3 C3 Shades Series	Keyed to Vita [®] 3D Master [®] Shades ²
A4	A4 B4	C3-418	4L2.5
		C3-420	4M2
		C3-421	4M3
		C3-423	4R2.5
		C3-525	5M2
B2	B1, B2	C3-526	5M3
		C3-203	2L1.5
C1	C1, C2 D2	C3-205	2M1
		C3-208	2R1.5
		C3-312	3M1
C3	C3, C4	C3-310	3L1.5
		C3-417	4L1.5
		C3-524	5M1

1. Vita[®] and 3D-Master[®] are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co. 2. Data on file.

1. Be sure the die is clean
2. Reapply the die release and allow to dry
3. Place the coping on the die and apply the Final Margin porcelain
4. Allow to dry several minutes, then lift the coping off the die. DO NOT ALLOW TO DRY COMPLETELY before lifting off the die.
5. If using an overglaze, carefully apply the glaze to the restoration.
6. Fire the restoration according to the recommended temperatures.

DENTSPLY

CERAMCO[®]3

Stains

GENERAL PRINCIPLES

1. Mix the DENTSPLY Stains and Overglaze with the DENTSPLY Stain and Glaze Liquid.
2. Apply the stain to the desired area.
3. Fire the porcelain on honeycomb sagger trays with metal pins.
4. Fire according to recommended temperatures.

TINT STAINS

The tint stains are designed to closely approximate the average range of tooth colors found in natural dentition and commercially available shade guides. Predetermined, ready to use tint stains are designed to increase chroma one or more shades in a given range, or subtly characterize and intensify neck and interproximal areas.

SHADE RANGES

i1 - Orange Tint stain	A1 - Orange Tint stain	C1 - Brown Tint stain
i3 - Orange Tint stain	A2 - Orange Tint stain	C2 - Brown Tint stain
i5 - Orange Tint stain	A3 - Orange Tint stain	C3 - Brown Tint stain
i7 - Orange Tint stain	A3.5 - Orange Tint stain	C4 - Brown Tint stain
i2 - Yellow Tint stain	A4 - Orange Tint stain	
i4 - Yellow Tint stain	B1 - Yellow Tint stain	D2 - Orange Tint stain
i6 - Yellow Tint stain	B2 - Yellow Tint stain	D3 - Yellow Tint stain
i8 - Yellow Tint stain	B3 - Yellow Tint stain	D4 - Orange Tint stain
	B4 - Yellow Tint stain	

APPLYING CERAMCO3 STAINS

1. Clean the porcelain surface by lightly blasting with 50 micron aluminum oxide at 20 psi, followed by cleaning with a steam cleaner or in a ultrasonic cleaner with distilled water for 5 minutes.
2. Use the tip of the staining brush and spread a thin film of stain liquid and overglaze mixture on the surface to be stained. If you wish to use a natural glaze, use only the glaze medium. This is done to simulate a glaze and enable the ceramist to judge the necessary correction.
3. Mix the stains with the stain medium to a thin consistency and apply the mixture to the porcelain surface. Do not allow the mixture to puddle. Apply the stain until the desired effect is achieved.

4. Fire according to recommended temperatures.

Applying the Stain

5. If stains are not required, the restoration can either be natural glazed, overglazed or mechanically polished using silicone wheels, porcelain polishing wheels, and flour pumice or aluminum oxide polishing paste. Diamond polishing paste may also be used but is not required.

Basic Build-Up Technique

Completed Illuminé Shaded Crown

Completed Crown with Internal Characterization

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Coloraciones		Pasta base	Cristales que se usan después de la pasta base	Cristales que se usan después de la pasta coloreada
i1, i2, i3, i4, i5, i6, i7, i8		Illuminé White	White	White
A1, A2, A3 A3.5, A4		Regular Regular	Medium Dark	Light Light
B1, B2 B3, B4		Regular Regular	Medium Dark	Light Light
C1 C2, C3, C4		Regular Regular	Medium Dark	Light Light
D2, D4 D3		Regular Regular	Medium Dark	Light Light
Ceramco®3 La serie de colores C3	Codificado en Vita® 3D-Master® Shade Series®2	Pasta base	Cristales que se usan después de la pasta base	Cristales que se usan después de la pasta coloreada
C3-101	1M1	Illuminé White	Medium	Light
C3-102	1M2	Illuminé White	Medium	Light
C3-203	2L1.5	Regular	Medium	Light
C3-204	2L2.5	Regular	Medium	Light
C3-205	2M1	Regular	Medium	Light
C3-206	2M2	Regular	Medium	Light
C3-207	2M3	Regular	Medium	Light
C3-208	2R1.5	Regular	Medium	Light
C3-209	2R2.5	Regular	Medium	Light
C3-310	3L1.5	Regular	Dark	Light
C3-311	3L2.5	Regular	Dark	Light
C3-312	3M1	Regular	Medium	Light
C3-313	3M2	Regular	Medium	Light
C3-314	3M3	Regular	Dark	Light
C3-315	3R1.5	Regular	Dark	Light
C3-316	3R2.5	Regular	Dark	Light
C3-417	4L1.5	Regular	Dark	Light
C3-418	4L2.5	Regular	Dark	Light
C3-419	4M1	Regular	Dark	Light
C3-420	4M2	Regular	Dark	Light
C3-421	4M3	Regular	Dark	Light
C3-422	4R1.5	Regular	Dark	Light
C3-423	4R2.5	Regular	Dark	Light
C3-524	5M1	Regular	Dark	Light
C3-525	5M2	Regular	Dark	Light
C3-526	5M3	Regular	Dark	Light

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